

PH-8900

ALMIGHTY OXYGEN JET



Please read through the operation manual before using the machine !

CONTENTS

USING STATEMEN	03
PRODUCT SUMMARIZE	03
IMPORTANCE OF OXYGEN TO HUMAN	04
WORKING THEORY	04
DISEASE SCHEME BY ANOXIC CELL	05
ACCESSORIES AND INSTRUCTION	06
MAIN FUNCTIONS	06
TO WHOM IS SUITABLE	07
TESTING	07
STEPS OF GENERAL CARE	08
IMPORTANCE OF FACIAL OXYGEN TREATMENT	08
STEPS OF OXYGEN SPRAY, OXYGEN INJECTION AND BREATH	09
SECTIONS OF FACIAL OXYGEN CARE AND OPERATIONS	10
ASSISTANT TREATMENTS	10
MAIN CHARACTERS	11
6-WEEKS TREATMENT EFFECTS SCHEME	11
CAUTIONS DURING OPERATION	11
MAINTENANCE AFTER USE	12
AFTER SERVICE AND GUARANTEE	12
THE CLIENT AS FOLLOWING IS NOT SUITABLE FOR THE TREATMENT.....	13
PRODUCT SERVICES AND MAINTENANCE	14

USING STATEMENT

A safe operation is contained in the manual. The information and pictures are the latest when the manual is published. The right of revision is reserved by the factory.

Caution: Never refit any part of the machine or use with products not belong to us, otherwise we will not offer the guarantee and be irresponsible for the loss caused by unallowed operations or actions.

- Please read through the instruction before using the machine!
- The instruction is a part of the product.
- Keep the instruction well for reading in future.

PRODUCT SUMMARIZE

The machine adopts the advanced transformation adsorption theory and separates the pure axenic oxygen with high density from the air under normal temperature. It solves the anoxic problem completely by oxygen breath, oxygen injection and atomizes oxygen spray. The oxygen injection treatment uses the pure oxygen to inject active product into the skin's bottom through injection tools with certain pressure and transfer pure oxygen into the deep-seated skin by pulse drive method.

Oxygen treatment is regarded as a versatile beauty treatment technique used widely in recent 10 years. Not only in iatrology history but also in beauty history, oxygen treatment made great contribution and was popular in all circles. In the field of iatrology research, the treatment has reached the painless state beyond the limit of micro-wound. The painless technique injected the bacterin or insulin into the body and absorbed by blood vessel. However, it's not an applied method for skin care or beauty using.

Pure Oxygen Injection System provides the painless treatment which makes use of oxygen to inject special active elements deep into the skin, such as oxygen serum and so on. These special active elements can reach the interspace between the scarfskin cells and make positive effects.

Pure Oxygen Injection System is able to inject 98% pure oxygen and active nutrition into the skin bottom by 2BAR pressure, which can accelerate the cell's metabolism and smoothen the skin to be in the bloom of youth.

IMPORTANCE OF OXYGEN TO HUMAN

In ancient ages, the oxygen content in the air is 38% to 50%. At the beginning of 18 century; the percentage fell to 30%. In the 20 century, the rate even fell to 23%. Up to now, the rate is less than 20%, which broke the minimum.

Healthy skin needs enough nutrition, but the air pollution gets more and more serious. As the oxygen content in the atmosphere is falling constantly, the lack of oxygen in the skin leads a lot of problem: oxygen lack in the skin, metabolism slows down, sensitive skin, sallow complexion, black spot, wrinkles and skin aging. Meanwhile, it will destroy the balance of lipid excreting and darken the skin. The research points that active elements and oxygen's movement can go beyond the interspace between cells and refresh the outer circulatory system for rebirth skin.


Oxygen and human's body are closely linked together. What caused our sickish body is the lack of oxygen in atmosphere, busy city life, long time facing to the computer, dry air office, lack of aerobics and bad eating habits. The researches show that we should breathe at least 1.44 gallon oxygen to keep the growth and rebirth of the 60 million cells in human's body.

WORKING THEORY

Using high pressure separate technique, it can make 98% pure oxygen from the air by ionization. Cooperate with special designed spray orifice, spray gun and mask, pure oxygen can be injected into the creasy and dry parts of the skin by oxygen injection, oxygen spray and oxygen smelling. In this way, it can make skin absorb oxygen fleetly; improve the cell's activity and metabolism and other unexpected effects.

DISEASE SCHEME BY ANOXIC CELL

※ We sump up the diseases as following



ACCESSORIES AND INSTRUCTION

1. Accessories



A. nutrition gun B. injection gun C. breath mask

2. Instruction

A. nutrition gun -----press the controller to release oxygen, pull back the controller to release distillate, atomized distillate and oxygen is available. There's a packet for distillate on gun head.

B. injection gun ----- used during oxygen injection

C. breath mask ----- used for oxygen breath and aroma treatment

MAIN FUNCTIONS

1. Absorb treatment

With concentrated activated distillate, pure oxygen massage can cure various pruritus, skin allergy, neurodermatitis, allergic dermatitis, eczema, sore skin and sunlight symptom

2. Injection treatment

By hypospray injection, it can inject the pure oxygen and beauty products into the cuticular bottom to intensify cell and its resistibility. What's more, it can expel rubbish easily; improve metabolism and many skin problems caused by lack of oxygen.

3. Smelling treatment

■Put on the mask and make effects with 98% oxygen and aromatic distillate as following

A. Breath system——nasosinusitis, faucitis, cold, etc;

B. Nervous system——neurasthenic, insomnia, migraine, etc;

C. Blood——leucopenia;

D. Digestive system——constipation, stomach disorders;

E. Metabolism & Incretion——diabetes, climacteric syndrome, menoxenia, dysmenorrheal, improve wound healing, metabolism and immunity, regulation of endocrine

F. Skin——various skin problems, eventual senility, improve blood circulation and metabolism, frap skin, dispel wrinkle and splash, whiten skin, balance lipin's excreting and improve sensitive skin

▲Have a rest after 15 minutes oxygen smelling.

TO WHOM IS SUITABLE

Pure Oxygen System is especially suitable for smoking women or women with problem skin. It can improve blood circulation and metabolism; improve the sensitive skin to the best state.

You will feel healthy and easy of your skin. The effect will be better if you treat with our company's distillate.

TESTING

High pressure oxygen with active elements is selected by dermatology research as experiment sample. The results are:

- ◆Oxygen injection is no harm to human's skin;
- ◆Active oxygen can reach deeply into the skin;
- ◆Active oxygen can quickly accumulates at the bottom of skin;
- ◆It can provide active oxygen and nutrient to the skin;
- ◆The treatment is operated by numerical control, which's safer than manual control.

STEPS OF GENERAL CARE

※ FACE

- 1、Cleaning: Wait for skin's dryness after cleaning; (Better with dermabrasion)
- 2、Oxygen spray: Spray oxygen on the face with "nutrition gun" to make the skin be used to the oxygen environment and improve its active cells and metabolism. (10 minutes)
- 3、Oxygen injection: Daub distillate on the face and use spray gun to inject oxygen by method of massage, lifting and skin scrubbing. (30 minutes)
- 4、Oxygen suction: Wear the oxygen mask and do oxygen suction care for about 15 minutes;
- 5、Surface application: Just use some moisture nutrition mask on the face;
- 6、Period of treatment: 12 times/1 course, twice every week in first 8 times, once every week in left 4 times.

IMPORTANCE OF FACIAL OXYGEN TREATMENT

Oxygen treatment is the healthiest beauty care. There's no blood vessel on the scarfskin. All the oxygen and nutrient should be transferred by the microcirculation system of the corium fat, which gets less and less when a person grows. Scientist found that the oxygen in the scarfskin fat when one is middle-aged is half of when one is young. If we can inject the oxygen into the skin's bottom, the cell will be activated and metabolism be improved. Oxygen treatment will make skin younger and easier.

STEPS OF OXYGEN SPRAY, OXYGEN INJECTION AND BREATH

1、Oxygen spray (Use nutrition gun)

Put some beauty liquid into the nutrition gun (according to customer's skin type), push back to spray nutrition and push down to spray oxygen.

ATTENTION: During oxygen spray on face, the pressure should be 5 LPM.

2、Oxygen injection (Use injection gun)

Choose right distillate to right skin. Dry skin can be moistened and whitened; Oily skin's pores can be shrunk and oil can be controlled well.

ATTENTION: Turn the pressure to 2LPM during oxygen injection. The injection time is 30 minutes for whole face from the forehead to chin. Please never inject on eyes. There are heads of 3 sizes of the injection gun.

3、Oxygen suction (Use breath mask)

Turn on the oxygen output function and wear the mask. Do surface application after oxygen suction.


※How will human be if there's no oxygen?

Brain anoxic——dizzy, headache, tinnitus, naupathia, nervous, breathe hard, blurry consciousness to death;

Skin anoxic——weaken update ability of skin, coarseness, largen pore, wrinkles, splash, black rim of eye, sallow complexion, etc.

ATTENTION: Turn the pressure to 3LPM during oxygen suction.

SECTIONS OF FACIAL OXYGEN CARE AND OPERATIONS



Make skin care according to the arrowhead and order above !


ASSISTANT TREATMENTS

It will be better when you use the oxygen injection system with assistant treatment such as soft light, color light, laser, micro-dermabration and other care

MAIN CHARACTERS

- It gathers oxygen atomization, oxygen injection and suction;
- It can provide 98% pure oxygen continuously;
- It makes the skin absorb nutrient easily;
- You will see unexpected effects if you use the machine with our special distillate;
- Nutrition gun is made of stainless steel and has a packet for distillate on gun head.

6-WEEKS TREATMENT EFFECTS SCHEME



An obvious result can be shown from above scheme that the skin becomes slicker and tender after 6 weeks' treatment.

CAUTIONS DURING OPERATION

- 1、 Attention must be paid to the operation condition:
 - (1) Never use the machine in dank or high temperature places to avoid creepage and operation malfunctions
 - (2) Operation condition: power 220 v .
- 2、 Attention must be paid before operation:
 - (1) Connect the power supply
 - (2) Confirm all the connections have been done completely
 - (3) Confirm if all the buttons are available by pressing a few times
 - (4) Confirm if all the functions and fittings are working correctly
 - (5) Confirm if the machine can start normally after no use for days

MAINTENANCE AFTER USE

Pure oxygen injection apparatus is an exact instruments which is not suitable in moist place. So, please choose the right storage place for it, especially for the power line.

- 1、 Cut off the power supply after use and pull out the pin;
- 2、 Avoid the sunlight's straight shining and moist places;
- 3、 Clean the fittings and put back after operation;
- 4、 Keep the machine and power line in different places separately after use. Especially, the power line should be kept in drafty place;
- 5、 Never place the used power line on top of the machine ;
- 6、 Avoid operating or keeping the machine in powerful humid place, such as bathroom ;
- 7、 Avoid short circuit and electric shock due to getting wet.

AFTER SERVICE AND GUARANTEE

- Our product has been inspected strictly before delivery. If you find any problems during operation, read the cautions again and repeat operation;
- If you meet problems or incomprehension, you should contact our agent in your place instead of fixing by yourself;
- Never split, repair or reform the machine privately except the manufacturer;
- Note the malfunction and details when you submit a repair.
- Guarantee regulations :

A、 One-year free guarantee of the mainframe is offered from the date you bought our product, except the man-caused malfunctions. (Fittings are not included)

B、 Attention! The followings are man-caused malfunctions:

- a: malfunction caused by incorrect operation;
- b: malfunction caused by incorrect reform;
- c: malfunction caused by natural disasters, such as flood, earthquake;
- d: malfunction caused by other situations.

THE CLIENT AS FOLLOWING IS NOT SUITABLE FOR THE TREATMENT.

CAUTIONS:

- 1, People who have heart disease, cardiac rhythm device configuration, high blood pressure.
- 2, Patients who has acute inflammation, deep vein thrombosis, thyroid, asthma, cancer and other .
- 3, Vascular bleeding, and trauma, the wound bleeding people.
- 4, Symptoms of skin inflammation, edema of the human body.
- 5, During pregnancy or within 3 months postpartum, it is not recommended for treatment of this equipment.
- 6, The body treatment is required 1 hour after meals to avoid feeling discomfort.
- 7, Apparatus should be thoroughly cleaned after use, so that equipment will last long, It increases the life span of equipment.
- 8, Do not unauthorizedly demolish, repair, alter the equipment (This will lead the failure formation or the accident).
- 9, At the end of using equipment, be sure to turn off the power switch and unplug the power plug.
- 10, Pay attention to fire prevention, protection, moisture-proof, anti-collision of the product, accessories should be properly placed.